

Portable Oxygen Concentrator 2-6L

This portable oxygen concentrator is designed for both travel and home use, providing supplemental oxygen. It features a split structure for easy installation and offers an oxygen flow rate of 2-6L.



ADDITIONAL IMAGES

PULSE OXYGEN SUPPLY
5 Levels Of Oxygen Supply Can Be Switched Freely

Free Switch

TRAVEL & HOME OXYGEN CONCENTRATOR
The Oxygen Concentration Is 90% ± 3%

- Split structure, easy installation
- Oxygen concentration 90% ± 3%
- On board 2L high flow
- Pulse oxygen supply
- Can be customized
- Current overload protection
- Charge by day and hour
- Complimentary Backpack
- Atomization function

Product Overview

HOME OXYGEN CONCENTRATOR
The Oxygen Concentration Is 90% ± 3%

- Split structure, easy installation
- Oxygen concentration 90% ± 3%
- On board 2L high flow
- Pulse oxygen supply
- Can be customized
- Current overload protection
- Charge by day and hour
- Complimentary Backpack
- Atomization function

APPLICATION SCENARIO
Widely Used

THE ELDERLY PREGNANT WOMEN
MOUNTAINEERING DRIVERS

Versatile Portable Oxygen Solutions

This portable oxygen concentrator provides a flexible solution for both home and travel use, offering an adjustable flow rate from 2L to 6L per minute. Engineered with advanced molecular sieve technology, it ensures high-purity oxygen delivery while maintaining energy efficiency. The device features a user-friendly interface with real-time monitoring of oxygen concentration and usage time, making it suitable for a wide range of users including the elderly and those requiring supplemental oxygen on the go.

Key Performance Metrics

Core Specifications

6 L/min Max Flow Rate	90 % Max Oxygen Concentration	120 W Rated Power Consumption	42 dB Noise Level
---------------------------------	---	---	-----------------------------

Technical Specifications

Output Pressure	0.03 - 0.07 Mpa
Packing Size	415 x 396 x 235 mm

Key Features

2-6L ADJUSTABLE FLOW

Oxygen concentration 90% ± 3%



Oxygen Concentration by Flow Rate

Flow Rate	Oxygen Concentration
2L	93%
3L	73%
4L	63%
5L	53%
6L	38%

Highlights

Pulse Oxygen Supply, Negative Ion Generator, Oil-free Compressor, SOS One-touch Button, Real-time Monitoring, Remote Control Included

Safety & Compliance

CARE AND SAFETY

To Be Safe When You Use Your Therapy Devices
Follow These Tips

- 1. Don't smoke. Tell others not to smoke around you, and stay away from smoky spaces when you go out. More oxygen in the air can make it easier for things to catch on fire. This also means that oily candles and matches can be dangerous.
- 2. Keep oxygen devices at least 3 feet away from stoves, cooktops, open flames, heaters, and electronic devices.
- 3. Avoid flammable liquids or cleaners when you use your device.
- 4. Use water-based moisturizer, not oil-based, on your lips, mouth, or nostrils.
- 5. Don't drink alcohol or take anything that might affect your breathing rate while using oxygen.
- 6. Talk with your doctor or nurse about how to clean your cannula, face mask, humidifier bottle, or catheter. Clean and replace them as recommended by your doctor.
- 7. Don't change your device's flow rate unless your doctor tells you it's OK.

Built-in Safety Protections

- Current overload protection
- Compressor high-temperature automatic shutdown
- Low oxygen concentration alarm
- Power failure alert